

Diaclean

Diaclean is a cleaning agent intended to clean the counting aperture of the DL3 and DLB5 haematology analysers

REF	Content
H19170	- 250 mL sodium hypochlorite solution

For in vitro use only.

GENERAL INFORMATION

Method	Manual cleaning
Shelf life	3 years from date of production
Storage	+15-25 °C

INTENDED USE

Diaclean is a cleaning agent intended to clean the counting aperture of the DL3 and DLB5 haematology analysers.

DIAGNOSTIC SIGNIFICANCE

For accurate measurements, it is essential to keep the aperture clean, therefore using the Diaclean solution is recommended after every 1000 cycles.

REAGENT COMPOSITION

Sodium hypochlorite (<5%)
 Distilled water

REAGENT PREPARATION

Diaclean is ready to use.

STORAGE AND STABILITY

- Diaclean:
- must never be subjected to temperatures below 15°C and above 25°C.
 - can be stored for 120 days (16 weeks) after opening.

- can be stored for 3 years (see expiration date on the product label) since the date of manufacture.

WARNINGS AND PRECAUTIONS

- For in vitro use only.
- Please read the Safety Data Sheet (MSDS) before use.
- Do not use directly after transport or directly after handling. Wait a few minutes.
- Cleaning solution can cause irritation of the eyes, skin and mucous membranes.
- In case of contact, rinse with plenty of water and seek immediate medical attention.
- Do not use the cleaning solution beyond the expiry date or beyond the open bottle time.
- Do not transfer / mix cleaning solutions whether they are of the same type or batch.
- Do not reuse an empty container as this may distort the results or damage the machine.
- Do not use the product in case of deterioration of the protective packaging.
- If the cleaning solution has been accidentally spilled, wipe it off immediately with a damp cloth.
- Do not use the product if there is any sign of deterioration (turbidity, color change, etc.)

TEST PROCEDURE

A detailed description is available in the user manual of the relevant analyser.

WASTE MANAGEMENT

Waste is disposed of according to local regulations.

USED SYMBOLS



Dispose of the tests and packaging appropriately

