


## DIAQUICK TROPONIN T

### GENERAL INFORMATION

- Cardiac Troponin T (cTnT), currently recognized as the most valuable diagnostic index for myocardial injury, has shown broad application prospects and has replaced creatine phosphate kinase MB isoenzyme (CK-MB) as the "gold standard" for the assessment of myocardial injury, especially for diagnosis of acute myocardial infarction.
- It plays an important role in the diagnosis of heart failure, unstable angina pectoris, myocarditis, drug-induced myocardial injury, cardiac injury monitoring in thoracic surgery, various critical diseases and multiple organ failure.
- The DIAQUICK Troponin T Cassette (Whole Blood/Serum/Plasma) is a simple test that utilizes a combination of anti-cTnT antibody coated particles and capture reagent to detect cTnT in whole blood, serum or plasma.
- The minimum detection level is 0.5 ng/mL.

### TEST PROCEDURE



### PRODUCT INFORMATION

Sample material	WB, Serum, Plasma
Shelf life	24 months
Method	Chromatographic immunoassay
Storage	2-30°C
Results	After 10 minutes



### ORDER INFORMATION

REF	Name	Specimen Type	Tests/Kit
Z20300CE	DIAQUICK Troponin T Cassette	WB, Serum, Plasma	30