

# DIAQUICK COVID-19 IgG/IgM

DIAQUICK COVID-19 CASSETTE

## GENERAL INFORMATION


- The 2019 new coronavirus pneumonia (NCP), or “COVID-19”, was discovered for the 2019 Wuhan Viral Pneumonia case in China and was named by the World Health Organization on January 12, 2020.
- Coronaviruses are enveloped RNA viruses that are distributed broadly among humans, other mammals and birds and cause respiratory, enteric, hepatic and neurological diseases.

## DIAQUICK COVID-19 IgG/IgM FEATURES

- A solid phase immunochromatographic assay
- A rapid, qualitative and differential detection of IgG and IgM antibodies to 2019 Novel Coronavirus
- Can be performed with serum, plasma and whole blood
- Fast results after 2-10 minutes
- Simple, time-saving procedure
- Little specimens, only 5µL of serum/plasma or 10µL of whole blood
- High sensitivity and specificity



## TEST PROCEDURE AND INTERPRETATION



## ORDER INFORMATION

REF	Name	Specimen Type	Tests/Kit
Z20101CE	DIAQUICK COVID-19 IgG/IgM Cassette	whole blood, serum, plasma	25