

**STATEMENT**  
**TO WHOM IT MAY CONCERN**

We



Dialab Produktion und Vertrieb von chemisch-technischen Produkten und  
Laborinstrumenten Gesellschaft m. b. H.  
IZ-NOE Sued, Hondastrasse, Objekt M55  
A-2351 Wiener Neudorf, Austria

Declare that:

DIAQUICK COVID-19 Ag Cassette is designated to detect nucleocapsid protein (NP) of SARS-CoV-2 virus. Epitope of capture antibody and detection antibody were not in the mutation area. Thus, above mentioned spike protein mutation SARS-CoV-2 and nucleocapsid protein mutation SARS-CoV-2 of Alpha(B.1.1.7), Beta(B.1.351), Gamma(P.1), Delta(B.1.617.2), Lambda(C.37), Zeta(P.2), Kappa(B.1.617.1), Eta(B.1.525), Iota(B.1.526), Epsilon(B.1.427), Theta(P.3) and Omicron(B.1.1.529) are not affecting virus performance of DIAQUICK COVID-19 Ag Cassette.



**DIALAB**  
Produktion und Vertrieb von chemisch-technischen  
Produkten und Laborinstrumenten Gesellschaft m.b.H.  
A - 2351 Wr. Neudorf, Z.N.O. Sued, Hondastr. Obj. M55  
Phone: +43(0)2236 660910 | Fax: +43(0)2236 660910 - 30  
E-Mail: office@dialab.at | Website: www.dialab.at

Wiener Neudorf, 2021-12-03

*2021-12-03*  
Vanessa Koch  
Head of Product Management