



BLOOD GROUPING & RHESUS TYPING



Features and Benefits

- **Agglutination method**
- **High specificity, high avidity and high potency**
- **Tested against NIBSC Reference Standards**
- **No instruments or extensive equipment required**
- **Manual test procedure**
- **Visual result interpretation**

ABO Grouping

Product	Clone / Cell line	Content	REF	CE
Anti-A, monoclonal	9113D10	10 mL 1000 mL	B05405 B08405	✓
Anti-B, monoclonal	9621A8	10 mL 1000 mL	B05406 B08406	✓
Anti-AB, monoclonal	152D12 + 9113D10	10 mL 1000 mL	B05407 B08407	✓
Anti-A1 Lectin	Dolichos Biflorus	5 mL	B05177	

Rhesus Typing

Product	Clone / Cell line	Content	REF	CE
Anti-D (IgG/IgM), monoclonal	RUM-1 + MS-26	10 mL 1000 mL	B05408 B08408	✓
Anti-C, monoclonal	MS-24	5 mL	B06411	✓
Anti-E, monoclonal	MS-258	5 mL	B06412	✓
Anti-c, monoclonal	MS-33	5 mL	B06413	✓
Anti-e, monoclonal	MS-16 + MS-63	5 mL	B06414	✓
Anti-Cw, monoclonal/	MS-110	2 mL	B13001	✓



Rare Blood Grouping Reagents

Product	Clone / Cell line	Content	REF	CE
Anti-H Lectin	Ulrex europaeus	2 mL	B09934	✓
Anti-K, monoclonal	MS-56	10 mL	B09909	
Anti-k, monoclonal	Human serum	2 mL	B18908	✓


Auxiliary Reagents

Product	Clone / Cell line / Description	Content	REF	CE
Anti-HG, polyspecific/monoclonal	Anti-IgG, BRIC-8	10 mL	B05181	✓
BSA 22%	Bovine serum albumin	10 mL	B05176	✓
BSA 30%	Bovine serum albumin	10 mL	B05175	✓
LISS Solution	Low ionic strength solution	10 mL 250 mL	B05178 B12250	✓
Papain	Stabilised papain	10 mL 1000 mL	B08000 B08001B	✓
Anti-D Negative Control	Anti-D reagent without antibodies	10 mL	B09936	✓

OEM / Bulk

Many of the listed products in our catalogue are also available in bulk, unlabelled and OEM version. For special requests or further information please contact us.





Austria (Head Office)

DIALAB® GmbH
IZ NOE-Sued
Hondastrasse, Objekt M55
A-2351 Wiener Neudorf
Phone: +43-2236-66 09 10-0
Fax: +43-2236-66 09 10-30
E-mail: office@dialab.at
www.dialab.at

China

DIALAB® ZGJ.
Building A, Room 211
No. 36, Huada Road, Zhangjiagang, Suzhou
Jiangsu Province, China
Phone: +86-0512-58711903
E-mail: office@dialab.at
www.dialab.at

Czech Republic

DIALAB® spol. s r.o.
Náměstí Osvoboditelů 1
CZ-153 00 Praha 5 - Radotín
Phone/Fax: +420-2-579 10 255
E-mail: office@dialab.cz
www.dialab.cz

Hungary

DIALAB® Kft.
Pasaréti út 31
H-1026 Budapest
Phone: +36-1-212 25 05
Fax: +36-1-335 73 70
E-mail: info@dialab.hu
www.dialab.hu

Serbia

DIALAB® d.o.o.
Carigradska 1
RS-11 008 Beograd 12
Phone: +381-11-32 38 745
Phone/Fax: +381-11-32 42 777
E-mail: dialab@sbb.rs

Slovak Republic

ALPHA DIALAB® s.r.o.
Sekurisova 16
SK-84102 Bratislava
Phone: +421-917 98 3711
E-mail: alpha@dialab.sk