

The logo for DIALAB, featuring the word "DIALAB" in a stylized, blue, outlined font with a registered trademark symbol (®) to the upper right.

DIALAB[®]

DIACHECK[®] PRO




Highlights

- Code free
- Wide test range
- Fast results
- For self-testing
- CE certified
- In compliance with EN ISO 15197:2015

General Features

| | |
|-------------------------|--------------------------------------|
| Methodology | Electrochemistry |
| Calibration | Plasma equivalent |
| Coding | Code free |
| Sample material | Capillary and venous whole blood |
| Test range | 10 – 900 mg/dL (0.6 – 50 mmol/L) |
| Hematocrit range | 20 – 60 % |
| Test time | 5 seconds |
| Sample volume | 0.5 µL |
| Memory | 1000 test results with date and time |



Order Information

| Product | Content | REF |
|---|---|--------|
| DiaCheck® Pro Blood Glucose Monitoring System Starter Package including Battery, Lancing Device and Carrying Case | 1 Blood Glucose Meter 10 Lancets 10 Test Strips | P13110 |
| DiaCheck® Pro Blood Glucose Test Strips | 2x 25 Test Strips | P13125 |
| DiaCheck® Pro Blood Glucose Control Solution | 1x 4 mL Negative Control 1x 4 mL Positive Control | P13134 |